

Technical Reference Guide

Determination of Protein Concentration

Introduction

Several methods are commonly used for determination of protein concentration. Measurement of the UV absorbance at 280 nm is most useful for pure protein solutions. Bradford and BCA assay methods are routinely used during protein purification and screening.

Protein Determination Using Absorbance at 280 nm

Determination of protein concentration by ultraviolet absorption (260 nm to 280 nm) depends on the presence of aromatic amino acids in proteins. Tyrosine and tryptophan absorb at approximately 280 nm. Higher orders of protein structure also may absorb UV light or modify the molar absorptivities of tyrosine and tryptophan, and thus the UV detection is highly sensitive to pH and ionic strength at which measurement is taken. Many other cellular components, and particularly nucleic acids, also absorb UV light. The ratio of A_{280}/A_{260} is often used as a criterion of the purity of protein or nucleic acid samples during their purification. The real advantages of this method of determining protein concentration are that the sample is not destroyed and that it is very rapid. Although different proteins will have different amino acid compositions and thus different molar absorptivities, this method can be very accurate when comparing different solutions of the same protein.

To make an accurate determination of protein concentration, you will have to produce a standard curve (A_{280}) with known amounts of purified protein. You will also have to provide a blank that is appropriate for the sample and contains the same concentrations of buffer and salts as the sample. It is often convenient to dialyze the sample and measure the absorbance of the sample (still in the dialysis sack) using the dialysis solution as the blank. Care must be taken to use quartz cuvettes, since glass absorbs UV light. A handy equation to estimate protein concentration that is often used is $[\text{Protein}] \text{ (mg/mL)} = 1.55 \cdot A_{280} - 0.76 \cdot A_{260}$

Protein Determination Using the Bradford Assay

The Bradford Assay, a colorimetric protein assay, is based on an absorbance shift in the dye Coomassie when bound to arginine and hydrophobic amino acid residues present in protein.

The (bound) form of the dye is blue and has an absorption spectrum maximum historically held to be at 595 nm. The anionic (unbound) forms are green and red. The increase of absorbance at 595 nm is proportional to the amount of bound dye, and thus to the amount (concentration) of protein present in the sample.

Unlike other protein assays, the Bradford Protein Assay is less susceptible to interference by various chemicals that may be present in protein samples.

Protein Determination Using the BCA Assay:

This assay measures the formation of Cu^{+1} from Cu^{+2} by the Biuret complex in alkaline solutions of protein. The first reaction occurs at lower temperatures and is the result of the interaction of copper and BCA with the following residues: cysteine, cystine, tryptophan and tyrosine. At elevated temperatures, it has been shown that the peptide bond itself is responsible for color development. This is why performing the assay at 37°C or 60°C increases the sensitivity and reduces the variation in the response of the assay to protein composition. The BCA reagent forms a complex with Cu^{+} , which has a strong absorbance at 562 nm.

Contact Information

North America

Customer Service: 800 638 8174 (toll free)
Scientific Support: 800 521 0390 (toll free)
scientific.support@lonza.com
Online Ordering: www.lonza.com

Europe

Customer Service: +32 87 321 611
Scientific Support: +49 221 99199 400
scientific.support.eu@lonza.com
Online Ordering: www.lonza.com

International

Contact your local Lonza Distributor
Customer Service: +1 301 898 7025, ext. 2322
Fax: +1 301 845 8291
scientific.support@lonza.com

International Offices

Australia	+61 3 9550 0883
Austria	0800 201 538 (toll free)
Belgium	+32 87 321 611
Brazil	+55 11 2069 8800
Denmark	+45 43 56 74 00
France	0800 91 19 81 (toll free)
Germany	0800 182 52 87 (toll free)
India	+91 22 4342 4000
Italy	+39 0363 45710
Japan	+81 3 5566 0612
Poland	+48 22 833 87 45
Singapore	+65 64914214
Spain	+34 902 531 366
Sweden	020 140 4410 (toll free)
Switzerland	0800 83 86 20 (toll free)
The Netherlands	0800 022 4525 (toll free)
United Kingdom	+44 118 979 5234

Lonza Cologne AG 50829 Cologne, Germany

The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. No statement is intended or should be construed as a recommendation to infringe any existing patent.

© Copyright 2009, Lonza Cologne AG. All rights reserved.
RG-DetProtConc 03/09
